

KERN OBB-A 1605

KERN

Infinity PH-Plan-objective 10 x



Category

Brand	Optics
Product category	Optical Instrument
Product group	PH-Objective

Construction

Anti fungus	✓
-------------	---

Environmental conditions

Storage temperature [Min]	-5 °C
Storage temperature [Max]	40 °C

Packing & Shipping

Readability force [d] (N)	1 d
Dimensions packaging (W×D×H)	65×40×40 mm
Net weight	0,3 kg
Shipping method	Parcel service
Net weight approx.	0,30 kg
Gross weight approx.	0,60 kg
Shipping weight	0,6 kg

Objective

Objective quality	Infinity Plan
Objective magnification	10 x
Numerical Aperture	0,25
Objective working distance	10,3 mm

Pictograms