

KERN OBB-A3203

KERN

Achromatic objective lens, 4 x /0,1 W.D. (27 mm)



Pictograms

STANDARD



Category

Brand	Optics
Product category	Optical Instrument
Product group	Microscope objective achromatic

Construction

Dimension (W×D×H)	27×25×25 mm
-------------------	-------------

Environmental conditions

Storage temperature [Min]	-5 °C
Storage temperature [Max]	40 °C

Packing & Shipping

Readability force [d] (N)	1 d
Dimensions packaging (W×D×H)	65×40×40 mm
Net weight	0,3 kg
Shipping method	Parcel service
Net weight approx.	0,30 kg
Gross weight approx.	0,60 kg
Shipping weight	0,6 kg

Objective

Objective quality	Achromatic
Objective magnification	4 x
Numerical Aperture	0,1
Objective working distance	27 mm